

## Simple Protocol for PSV-Based Neutralization Assay

### **SARS-CoV-2 (2019nCoV) pseudotype virus (pseudovirus, PSV) for COVID-19 related vaccines and neutralizing antibodies evaluation**

The outbreak of COVID-19, caused by SARS-CoV-2 (2019-nCoV), has been a global public health threat and caught the worldwide concern. Due to its high pathogenicity and infectivity<sup>1</sup>, live SARS-CoV-2 should be handled under biosafety level 3 (BSL-3) conditions. GeneMedi has developed SARS-CoV-2 pseudovirus production system, from which the SARS-CoV-2 pseudotyped virus can be handled in biosafety level 2 (BSL-2)<sup>2</sup>.

GeneMedi's SARS-CoV-2 (2019nCoV) pseudotype virus (pseudovirus, PSV) based neutralization assay is a standard evaluation procedure for COVID-19 related vaccines and neutralizing antibodies potency evaluation. GeneMedi's SARS-CoV-2 PSV is the core ingredient of diagnostics for neutralization serology after vaccinotherapy.

GeneMedi's SARS-CoV-2 pseudotyped virus includes wildtype and the spike mutation variants (D614G, S943P, V367F, G476S, V483A, H49Y, Q239K, A831V, P1263L, D839Y/N/E: D839Y, D839N, D839E). The GeneMedi's SARS-CoV2 PSV panel help for all-in-one vaccinotherapy evaluation.

### **Application**

#### **SARS-CoV-2(2019nCoV) Pseudotyped Virus Based Neutralization Assay<sup>3</sup>**

Coronavirus disease 2019 (COVID-19) pandemic is caused by SARS-CoV-2 (2019nCoV) infection, a newly emerged novel coronavirus spreading worldwide. Current efforts are focusing on development of specific antiviral drugs. Therapeutic neutralizing antibodies (NAbs) against SARS-CoV-2(2019-nCoV) will be greatly important therapeutic agents for the treatment of COVID-19. The availability of therapeutic NAbs against SARS-CoV-2 will offer benefits for the control of the current pandemic and the possible re-emergence of the virus in the future, and their development therefore remains a high priority.

GeneMedi's NAbs has been validated to reduce SARS-CoV-2 lentivirus-based pseudo virus infectivity and thereby blocking the entry of the Coronavirus to its effector/targeting cell: human ACE2-HEK293T cell (hACE2-HEK293T, Cat. GM-SC-293T-hACE201). GeneMedi's SARS-CoV-2 (2019nCoV) Nabs can act as a benchmark of neutralizing antibodies discovery against COVID-19.

GeneMedi's Pseudovirus Based Neutralization Assay (PBNA) is a conventional assay method that is suitable for High-Throughput Screening (HTS) without live virus engaged. The Pseudovirus Based Neutralization Assay can be used for evaluating

- 1) Neutralizing antibodies (NAbs)<sup>3,4</sup>
- 2) Peptides blockers<sup>5,6</sup> (peptide inhibitors) or protein<sup>7,8</sup>
- 3) Types of Vaccines (Immunized serum)<sup>9</sup>
- 4) Compounds targeting Spike induced cell-fusion<sup>10</sup>.

## Materials

1. SARS-CoV-2 Pseudovirus-RFP-fLuciferase ([GM-2019nCoV-PSV01](#))
2. Effector cell: Alternative
  - A. hACE2-HEK293T stable cell line ([GM-SC-293T-hACE2-01](#))
  - B. Wildtype HEK293T cell line, hACE2 vector for transfection ([GMV-V-2019nCoV-041](#))
3. Neutralizing antibodies (NAbs) ([GMP-V-2019nCoV-SnAb001~GMP-V-2019nCoV-SnAb005](#))

### **Pseudotyped virus of SARS-CoV-2 Spike Mutation Variants (D614G, S943P, V367F, G476S, V483A, H49Y, Q239K, A831V, P1263L, D839Y/N/E:D839Y, D839N, D839E)**

Catalog No.	Pseudotyped virus of SARS-CoV-2 Spike Mutation Variants
GM-2019nCoV-PSV02	<a href="#">Spike D614G mutation SARS-CoV-2(2019nCoV) Pseudotyped virus</a>
GM-2019nCoV-PSV03	<a href="#">Spike S943P mutation SARS-CoV-2(2019nCoV) Pseudotyped virus</a>
GM-2019nCoV-PSV04	<a href="#">Spike V367F mutation SARS-CoV-2(2019nCoV) Pseudotyped virus</a>
GM-2019nCoV-PSV05	<a href="#">Spike G476S mutation SARS-CoV-2(2019nCoV) Pseudotyped virus</a>
GM-2019nCoV-PSV06	<a href="#">Spike V483A mutation SARS-CoV-2(2019nCoV) Pseudotyped virus</a>
GM-2019nCoV-PSV07	<a href="#">Spike H49Y mutation SARS-CoV-2(2019nCoV) Pseudotyped virus</a>
GM-2019nCoV-PSV08	<a href="#">Spike Q239K mutation SARS-CoV-2(2019nCoV) Pseudotyped virus</a>
GM-2019nCoV-PSV09	<a href="#">Spike A831V mutation SARS-CoV-2(2019nCoV) Pseudotyped virus</a>
GM-2019nCoV-PSV10	<a href="#">Spike P1263L mutation SARS-CoV-2(2019nCoV) Pseudotyped virus</a>
GM-2019nCoV-PSV11	<a href="#">Spike D839Y/N/E-D839Y mutation SARS-CoV-2(2019nCoV) Pseudotyped virus</a>
GM-2019nCoV-PSV12	<a href="#">Spike D839Y/N/E-D839N mutation SARS-CoV-2(2019nCoV) Pseudotyped virus</a>
GM-2019nCoV-PSV13	<a href="#">Spike D839Y/N/E-D839E mutation SARS-CoV-2(2019nCoV) Pseudotyped virus</a>

## Protocol:

*If your effector cell is hACE2-HEK293T stable cell line, please begin in Step 2.*

1. Transfect HEK293T with hACE2-GFP vector ([GMV-V-2019nCoV-041](#)) 24hrs before planting the cell into 96-well.
2. Plant the hACE2-HEK293T into 96-well (5,000~10,000 per well) overnight before SARS-CoV-2 PSV infection.
3. Generation of 100ul PSV-Sample mixture:

100ul PSV-Sample mixture	Volume
GM-2019nCoV-PSV01*	50ul or 5ul
Sample(NAbs, peptides, serum, etc)	flexible (According to your own products)
Total	add culture medium to 100ul
* For GM-2019nCoV-PSV01-1, add 50ul in recommendation (range from 20ul~100ul). For GM-2019nCoV-PSV01-2, add 5ul in recommendation. (range from 2ul~10ul).	

Incubate PSV-Sample mixture for 1h at room temperature.

4. Remove the medium of effector cells in 96-well, add 100ul PSV-Sample mixture to 96-well for infection, 3 replicates per group.
5. Fluorescence imaging (RFP) 72hrs after SARS-CoV-2 PSV infection. The firefly luciferase reporter is measured following the Promega Luciferase Assay Reagent manual.

### Tips

#### If your samples are serum

A standard curve should be generated using serially diluted Nabs (neutralizing antibodies) as a positive control.

#### If your samples are therapeutic antibodies or peptides candidates

Dilute the samples into concentration gradient for IC50 value evaluation.

## Product information

<b>Catalot Number</b>	GM-2019nCoV-PSV01
<b>Products Name</b>	SARS-CoV-2 Pseudovirus-RFP-fLuciferase
<b>Reporter</b>	RFP+Firefly Luciferase
<b>Ligand-Receptor</b>	Spike-ACE2
<b>Effector cell recommended</b>	hACE2-HEK293T ( <a href="#">Cat.GM-SC-293T-hACE201</a> )
<b>Size</b>	0.5ml/vial
<b>Products description</b>	<p>GeneMedi's SARS-CoV-2 Pseudovirus-RFP-fLuciferase (GM-2019nCoV-PSV01) is recombinant pseudotyped lentiviral particles containing SARS-CoV-2 spike protein to mimic SARS-CoV-2 (2019nCoV) cell infection and cell entry.</p> <p>The SARS-CoV-2 pseudovirus particles encode firefly luciferase and RFP in their lentiviral vector genome. The firefly luciferase and RFP gene will be strongly expressed after the SARS-CoV-2 pseudovirus entry into ACE2-expressing cells. Actually, 293T-hACE2(human ACE2 overexpression stable HEK293T cell lines) is normally used as effector cell (GM-SC-293T-hACE2-01).</p> <p>GM-2019nCoV-PSV01 is a powerful tool for SARS-CoV-2 related vaccine efficacy evaluation, neutralizing antibodies, peptides blockers competitors neutralization assay, and tissue-specific infection determination.</p>
<b>Bioactivity validation</b>	Validated in hACE2-HEK293T Cell Entry
<b>Storage</b>	Store at -20°C to -80°C under sterile conditions. Avoid repeated freeze-thaw cycles.

# GeneMedi SARS-CoV-2 Pseudovirus (PSV) Based Cell Entry




Figure. GeneMedi-SARS-CoV-2 Pseudovirus (PSV) Based Cell Entry validation with fluorescence(A), FACS (B) and Luciferase activity (C) after 72 hours of HEK293T infection. hACE2 significantly enhanced the infection efficiency of the SARS-CoV-2 PSV.

GeneMedi SARS-CoV-2 PSV-Luciferase ([Cat.GM-2019nCoV-PSV01](#)) is recombinant pseudotyped lentiviral particles containing SARS-CoV-2 spike protein to mimic SARS-CoV-2 (2019nCoV) cell infection.

GM-2019nCoV-PSV01 is a powerful tool for SARS-CoV-2 related vaccine efficacy evaluation, neutralizing antibodies, peptides blockers competitors neutralization assay, and tissue-specific infection determination.

# GeneMedi COVID-19 neutralizing antibodies assay system

--Nab discovery and vaccines evaluation through SARS-CoV-2 wildtype/mutant variants pseudovirus based neutralizing assay(PBNA) and Spike-ACE2 competition binding assay

## GeneMedi-SARS-CoV-2 WT and Spike Mutation Variants Pseudovirus (PSV) Based Cell Entry





Figure. The Pseudovirus (PSV) Based Cell Entry assay was performed on 293T-hACE2 cells infected with [GeneMedi-SARS-CoV-2 WT and Spike Mutation Variants \(D614G, S943P, V367F, G476S, V483A, H49Y, Q239K, A831V, P1263L, D839Y/N/E:D839Y,D839N,D839E\) Pseudovirus \(PSV\)](#) Infection rate was determined by RFP+GFP+/GFP+ with FACS validation.

## GeneMedi's anti-2019-nCoV Spike Neutralizing antibodies (Nabs) and Spike RBD protein binding validation



**Figure.** The binding of GeneMedi's anti-2019-nCoV Spike Neutralizing antibodies (Nabs) to Recombinant 2019-nCoV(SARS-CoV-2) Spike RBD protein ([GMP-V-2019nCoV-SRBD001](#)) at 5.0ug/ml (100uL/well) was measured by ELISA.

- A. [GMP-V-2019nCoV-SnAb001](#): Anti-2019-nCoV Spike (Spike RBD domain) human monoclonal neutralizing antibody (IgG1)  
 B. [GMP-V-2019nCoV-SnAb002](#): Anti-2019-nCoV Spike (Spike RBD domain) human monoclonal neutralizing antibody (IgM)  
 C. [GMP-V-2019nCoV-SnAb003](#): Anti-2019-nCoV Spike (Spike RBD domain) human monoclonal neutralizing antibody (IgA)  
 D. [GMP-V-2019nCoV-SnAb004](#): Anti-2019-nCoV Spike (Spike RBD domain) mouse monoclonal neutralizing antibody (IgG1)  
 E. [GMP-V-2019nCoV-SnAb005](#): Anti-2019-nCoV Spike (Spike RBD domain) Cynomolgus monoclonal neutralizing antibody (IgG1)

## GeneMedi's anti-2019-nCoV Spike Neutralizing antibodies (Nabs) competitive binding assay validation




Figure. GeneMedi's anti-2019-nCoV Spike Neutralizing antibodies (Nabs) block Recombinant 2019-nCoV(SARS-CoV-2) Spike RBD protein ([GMP-V-2019nCoV-SRBD001](#)) and hACE2 ([GMP-H-ACE2002](#)) binding.

- A. [GMP-V-2019nCoV-SnAb001](#): Anti-2019-nCoV Spike (Spike RBD domain) human monoclonal neutralizing antibody (IgG1)
- B. [GMP-V-2019nCoV-SnAb002](#): Anti-2019-nCoV Spike (Spike RBD domain) human monoclonal neutralizing antibody (IgM)
- C. [GMP-V-2019nCoV-SnAb003](#): Anti-2019-nCoV Spike (Spike RBD domain) human monoclonal neutralizing antibody (IgA)
- D. [GMP-V-2019nCoV-SnAb004](#): Anti-2019-nCoV Spike (Spike RBD domain) mouse monoclonal neutralizing antibody (IgG1)
- E. [GMP-V-2019nCoV-SnAb005](#): Anti-2019-nCoV Spike (Spike RBD domain) Cynomolgus monoclonal neutralizing antibody (IgG1)



; YbYA YX]G5 FG!7 cJ!&K H'UbX'Gd]\_Y'Ai HU]cb'J U]Ubhg'DgYi Xcj ]fi g'fDGJ'L6 UgYX'BYi kfU]n]b[ ' 5 ggUmik ]A ; YbYA YX]fj'Ubh]!&\$% !b7 cJ'Gd]\_Y'BYi kfU]n]b[ 'Ubh]VcX]Yg'fBUVgk'




Figure. The Pseudovirus (PSV) Based Neutralizing Assay was performed on 293T-hACE2 cells infected with [GeneMedi-SARS-CoV-2 WT and Spike Mutation Variants \(D614G, S943P, V367F, G476S, V483A, H49Y, Q239K, A831V, P1263L, D839Y/N/E:D839Y,D839N,D839E\) Pseudovirus \(PSV\)](#) under treatment of GeneMedi's anti-2019-nCoV Spike Neutralizing antibodies (Nabs). Inhibition rate was determined by comparing the relative RFP+GFP+/GFP+ rate.

## References

- 1 XiaolongCai. An Insight of comparison between COVID-19 (2019-nCoV) and SARS-CoV in pathology and pathogenesis. *Preprint*, doi:10.31219/osf.io/hw34x (2020).
- 2 Jean K. Millet<sup>1</sup>, Tiffany Tang<sup>3</sup>, Lakshmi Nathan<sup>3</sup>, Javier A. Jaimes<sup>4</sup>, Hung-Lun Hsu<sup>3,5</sup>, & Susan Daniel<sup>3</sup>, G. R. W. Production of Pseudotyped Particles to Study Highly Pathogenic Coronaviruses in a Biosafety Level 2 Setting. *J Vis Exp*, doi:10.3791/59010 (2019).
- 3 Nie, J. *et al.* Establishment and validation of a pseudovirus neutralization assay for SARS-CoV-2. *Emerg Microbes Infect* **9**, 680–686, doi:10.1080/22221751.2020.1743767 (2020).
- 4 Jiang, S., Hillyer, C. & Du, L. Neutralizing Antibodies against SARS-CoV-2 and Other Human Coronaviruses. *Trends Immunol*, doi:10.1016/j.it.2020.03.007 (2020).
- 5 Zhang, G., Pomplun, S., Loftis, A. R., Loas, A. & Pentelute, B. L. The first-in-class peptide binder to the SARS-CoV-2 spike protein. doi:10.1101/2020.03.19.999318 (2020).
- 6 Xia, S. *et al.* Inhibition of SARS-CoV-2 (previously 2019-nCoV) infection by a highly potent pan-coronavirus fusion inhibitor targeting its spike protein that harbors a high capacity to mediate membrane fusion. *Cell Res* **30**, 343–355, doi:10.1038/s41422-020-0305-x (2020).
- 7 Monteil, V. *et al.* Inhibition of SARS-CoV-2 Infections in Engineered Human Tissues Using Clinical-Grade Soluble Human ACE2. *Cell*, doi:10.1016/j.cell.2020.04.004 (2020).
- 8 Lei, C. *et al.* Neutralization of SARS-CoV-2 spike pseudotyped virus by recombinant ACE2-Ig. *Nat Commun* **11**, 2070, doi:10.1038/s41467-020-16048-4 (2020).
- 9 HuajunBai, X., XiaolongCai. Vaccines And Advanced Vaccines: -A landscape for advanced vaccine technology against infectious disease ,COVID-19 and tumor. *Preprint*, doi:10.31219/osf.io/ypgx4 (2020).
- 10 Hoffmann, M. *et al.* The novel coronavirus 2019 (2019-nCoV) uses the SARS-coronavirus receptor ACE2 and the cellular protease TMPRSS2 for entry into target cells. *bioRxiv*, doi:10.1101/2020.01.31.929042 (2020).